


Anchorage process Quick Anchor

Anchor set on a ligament

REINFORCEMENT FIXATION SYSTEM



COUSIN
BIOTECH

Anchorage process Quick Anchor

Anchor set on a ligament Polyethylene terephthalate

DESCRIPTION :

The **anchorage process Quick Anchor** consists of an anchor \varnothing 3,2mm set on a polyester ligament L800mm with needles.

There are 2 types :

An anchor with a ligament \varnothing 0,6mm
FBIO AN306L



An anchor with a ligament \varnothing 0,8mm
FBIO AN308L



INDICATIONS :

The **anchorage process Quick Anchor** could be used in the surgery :

- for the shoulder,
- for the knee,
- and for the ankle.

as reinforcement fixation system.

CHARATERISTICS	FBIO AN306L	FBIO AN308L
Anchor	1 anchor \varnothing 3.2mm	1 anchor \varnothing 3.2mm
Braid	Polyester braid \varnothing 0.6mm L800mm	Polyester braid \varnothing 0.8mm L800mm
Needles	2 needles $\frac{1}{2}$ curve L30mm \varnothing 0.9mm	2 needles $\frac{3}{8}$ ^{ème} curve L50mm \varnothing 1.2mm
Breaking strength (daN)	8.5	27
Elongation (%)	11	17.5

PRESENTATION

Item numbers and sizes :



FBIO AN306L

FBIO AN308L

*Double packaging and rigid box. Sterilised products.
Manufacture and packaging made in controlled atmosphere, of class 10 US out of activity measured (norm FS 209E)*



DISTRIBUTED BY :

8, rue de l'abbé Bonpain – 59117 WERVICQ-SUD - FRANCE

Tel +33 (0) 3 20 14 40 00 – Fax +33 (0) 3 20 14 41 22

e-mail contact@cousin-biotech.com

www.cousin-biotech.com

